

C65T



Characteristics

Capacity	65 lt (2,30 cu ft)
External Dimensions (W x D x H)	cm 70 x 44 x 47
	<i>Internal Dimensions (W x D x H) cm 47 x 37,5 x 35 (+ additional compartment 15 x 37,5 x 17); Weight kg 24</i>
Structure and Insulation	Atoxic plastic material int./est., thickness 40 mm
Door	Nr 1 blind lid
Key Lock	
Internal equipment	Nr 1 open wire basket dim. cm. 46 x 36 x 32 (W x D x H) Nr.1 cable 12V
Controller	Digital electronic thermostat
Alarms	Acoustic and Visual for temperature deviations
Refrigeration Type	Static
Defrost	Manual

Technical data

Details	<ul style="list-style-type: none">• Refrigerant type: R134a CFC free• Mains Voltage 12V;• Nominal consumption 12 Vcc: 3,68 A• Nominal power consumption: 45 W• Ambient Temperature: +10/+32 °C; humidity 60%• Connectors: Socket for 12/24V
----------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

In accordance with	2014/30/CE, 2014/35/CE and following modifications
---------------------------	----------------------------------------------------

Accessories

Trasformer 230V	Setpoint +22°C; Settable between +30/-20°C
------------------------	--------------------------------------------

Accessories

Trasformer 230V

Graphic Chart Recorder

Datalogger

FIOCCHETTI SCIENTIFIC S.R.L.

Via Panagulis, 48 - 42045 Luzzara (RE) Italy
Tel. +39 0522 976232- Fax +39 0522 976028
info@fiocchetti.it - www.fiocchetti.it
P.I.00934960352
Registro A.E.E. IT8020000002833