

SUPERFREEZER ECO 70**Characteristics**

Temperature Range	-16/-32 °C
Capacity	70 lt (2,47 cu ft)
External Dimensions (W x D x H)	cm 55 X 66 X 63,6 <i>Weight kg 30</i>
Structure and Insulation	exterior sheet steel epoxy coated white colour/interior atoxic plastic material, thickness 65 mm
Door	Nr 1 blind door
Internal equipment	nr.2 drawers
Controller	Mechanical thermostat + LCD Digital Thermometer
Alarms	Visual for high temperature deviation
Refrigeration Type	Static
Defrost	Manual

Technical data

Details	<ul style="list-style-type: none">• Refrigerant type: R600a CFC free• Mains Voltage 230V/1/50Hz• Energy Efficiency Class: A+• Noise level < 40 dB• Ambient Temperature: +10/+43 °C; humidity 60%• Connector: Schuko
In accordance with	2014/30/CE, 2014/35/CE and following modifications
Package	Strong Cardboard on wooden pallet <i>Dim. cm.66x72x83; Kg.37</i> • Wooden crate on wooden pallet available: dim. cm.66x83x89; Kg.56

Accessories**Key Lock****DMLP UNIT**

External module with DMLP DIGITAL MONITOR: Monitor/alarm/recording system with rechargeable battery
- Temperatures and alarms recorder on removable Secure Digital (SD Back-Up memory): it can store 5 years data, recording temperature and alarms every minute. The SD card can be connected to any PC: a software (Capture) is included, so that you can download, save and print all data from your PC
- Battery backup: ensures alarms and recording of temperatures even in case of power failure. The battery system lasts about 12 hours.
- Dry contact for remote connection + RS485 port

DMLP PRINTER**FIOCCHETTI SCIENTIFIC S.R.L.**

Via Panagulis, 48 - 42045 Luzzara (RE) Italy
Tel. +39 0522 976232 - Fax +39 0522 976028
info@fiocchetti.it - www.fiocchetti.it
P.I.00934960352
Registro A.E.E. IT8020000002833