

SUPERFREEZER ECO 130



Characteristics

Temperature Range	-16/-28 °C
Capacity	130 lt (4,60 cu ft)
External Dimensions (W x D x H)	cm 55 X 66 X 86 <i>Weight kg 45</i>
Structure and Insulation	interior sheet steel epoxy coated white colour/interior atoxic plastic material, thickness 80 mm
Door	Nr 1 blind door
Internal equipment	nr.4 drawers
Controller	Mechanical thermostat + LCD Digital Thermometer
Alarms	Acoustic and Visual for temperature deviations
Refrigeration Type	Static
Defrost	Manual

Technical data

Details	<ul style="list-style-type: none">• Refrigerant type: R600a CFC free• Mains Voltage 230V/1/50Hz• Energy Efficiency Class: A++• Noise level < 40 dB• Ambient Temperature: +10/+32 °C; humidity 60%• Connector: Schuko
In accordance with	2014/30/CE, 2014/35/CE and following modifications
Package	Strong Cardboard on wooden pallet <i>Dim. cm.66x70x103; Kg.50</i> <ul style="list-style-type: none">• Wooden crate on wooden pallet available: dim. cm.75x76x106; Kg.71

Accessories

Key Lock

DMLP UNIT

External module with DMLP DIGITAL MONITOR: Monitor/alarm/recording system with rechargeable battery
- Temperatures and alarms recorder on removable Secure Digital (SD Back-Up memory): it can store 5 years data, recording temperature and alarms every minute. The SD card can be connected to any PC: a software (Capture) is included, so that you can download, save and print all data from your PC
- Battery backup: ensures alarms and recording of temperatures even in case of power failure. The battery system lasts about 12 hours.
- Dry contact for remote connection + RS485 port

DMLP PRINTER

FIOCCHETTI SCIENTIFIC S.R.L.

Via Panagulis, 48 - 42045 Luzzara (RE) Italy
Tel. +39 0522 976232- Fax +39 0522 976028
info@fiocchetti.it - www.fiocchetti.it
P.I.00934960352
Registro A.E.E. IT8020000002833