

B20W



Characteristics

External Dimensions (W x D x H) cm 39 x 52 x 64

Structure and Insulation stainless steel inside and outside

Kind of ice cubers

Motor cooling water

Production Kg/24h (up to) 24

Bin capacity (kg) 8

Weight kg 36

Technical data

Details

- Refrigerant type: R134a CFC free
- Mains Voltage: 230V/1/50Hz
- Power Consumption: 6 kW/h
- Water consumption l/h: 22
- Electrical consumption: 310 W
- Fuse (A): 10
- Cooling Capacity Btuh/W: 560
- Water temperature min/max: 5°C-35°C
- Pressure min/max: 1 bar/5 bar
- Noise level < 58 dB
- Ambient Temperature: +18/+43 °C; humidity 60%

In accordance with 2014/30/CE, 2014/35/CE and following modifications

Package Strong Cardboard on wooden pallet

Dim. cm.47 x 63 x 90; Kg.44

- *Wooden crate on wooden pallet available: dim. cm. 53 x 72 x 93; Kg. 68*

note Available also air cooled

Accessories

Kit 4 feet (+135mm H)

60 Hz Frequency

Wooden crate Package

FIOCCHETTI SCIENTIFIC S.R.L.

Via Panagulis, 48 - 42045 Luzzara (RE) Italy

Tel. +39 0522 976232- Fax +39 0522 976028

info@fiocchetti.it - www.fiocchetti.it

P.I.00934960352

Registro A.E.E. IT8020000002833